



TAG Sync x1 for BioExplorer


Theta-Alpha-Gamma Synchrony Single-Channel Solution

- 1 - Raw EEG for channel 1
- 2 - Spectral display for channel 1
- 3 - Amplitude vs time for Rewards 1 & 2
- 4 - Reward 1 amplitude (e.g., ultra-low frequency)
- 5 - Reward 2 amplitude (e.g., SMR)
- 6 - Multiple inhibit bins
- 7 - Status lights (e.g. phase reset candidates)



Key:

BioExplorer transport controls, menu and tools.



This is one of 10 screens in the 1-Channel TAG Sync Suite