



# TAG Sync x2 for Infiniti


Theta-Alpha-Gamma Synchrony

This is one of 15 screens included in this suite.



Key

- 1 - Raw EEG for Channel C
- 2 - Raw EEG for Channel D
- 3 - Spectral display for C+D
- 4 - Amplitude vs Time for Rewards 1 & 2
- 5 - TAG Nav can flag a candidate phase reset
- 6 - Typical Infiniti Lo & Hi inhibit instruments for C & D



- Reward 3 - Biofeedback (GSR, Temp, HRV)
- Reward 2 - EEG band such as alpha or gamma
- Reward 1 - EEG band such as delta or theta



Key

- 7 - Six programmable inhibit bins for Channel C
  - 8 - Six programmable inhibit bins for Channel D
- Note: % instrument below each bin indicates time above threshold.