


TAG Sync x2 for BioExplorer

Theta-Alpha-Gamma Synchrony

This is one of six screens in the TAG Sync 2 Suite

- 1 - Raw EEG for channel 1
- 2 - Raw EEG for channel 2
- 3 - Spectral display for 1+2
- Key: 4 - Amplitude vs time for Rewards 1 & 2
- 5 - Reward 1 thresholds
- 6 - Reward 2 thresholds

BioExplorer transport controls, menus and tools.



- 7 - Multiple inhibit bins for EEG channel 1
- 8 - Multiple inhibit bins for EEG channel 2
- Key: 9 - Status lights (phase reset candidates)